erbrich
hippo-dent®

Quality made in Germany

Pferdezahnschleifmaschine
Tooth Grinding Unit for Horses
Meuleuse pour dents de cheval
Fresa dentale
Torno dental

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65E-038.07
Reinigungsbürste, Messing


65E-500.28
Reinigungsbürste, Messing groß

65E**GB****Tooth Grinding Unit for Horses****erbrich
hippo-dent®****65E-169.23****Working attachment**

Reinforced, anodized aluminium shaft tube capable of withstanding extreme loads.
 Knurled section for a sure grip
 Rasp head made of high-alloy surgical stainless steel means no flaking chrome or nickel plating
 Small head for easier access to posterior molars
 Rasp head is protected by five hydrodynamic seals to prevent dirt from entering the gears
 New gearwheel geometry ensures efficient power transmission, quiet operation, low vibration, less friction and a long service life
 The gears are mounted in high-alloy bearings of the type used in the aircraft industry – a further guarantee of a long service life

Rasp disk

Main body made of hardened stainless steel for dimensional stability and to protect the gears
 Optimised geometric rasp disk holder for power transmission with maximum efficiency
 Diamond-coated as standard for long life - up to 300 treatment cycles
 Grinding with the diamond coating produces a texture on the tooth surface which is comparable to natural wear
 State-of-the-art coating technology prevents any flaking of the diamond coating

**Revolutionary power transmission system**

Newly designed connection between drive unit and attachment eliminates the need for a slip clutch – even simpler
 Direct drive via unbreakable shaft with multiple bearings, so there is no belt to snap or slip:
 no power loss, quiet operation, long service life

**Speed control**

12 speeds from 1100 to 2550 rpm (relative to the rasp disk) can be preselected and set for relaxed working
 The unit can be set to suit the individual tooth hardness, which means shorter treatment times

**Drive unit**

Low-maintenance 650 W professional power drill
 Drill stops automatically when carbon brushes are worn down
 Replaceable 6 m rubber-sheathed power cable for greater freedom of movement and no twisting

Technical data

230 V AC, 50 Hz
 Output 650 W
 Speed 0–2550 rpm (relative to the rasp disk)
 Doubly insulated motor – Class II
 Free from interference
 Conforms to CE marking requirements
 Weight 2.300 g

Scope of supply:

**65E-169.23, 65E-169.01, 65E-169.03,
 65E-038.07, 65E-038.08, 65E-038.11**

Accessories

65E-169.01
Spare diamond-coated disk
Long life (standard)

65E-169.02
Spare diamond-coated disk
High Speed Cutting

65E038.02
Spare diamond-coated disk
Long life (standard)
for EM2000

65E-038.09
Spare diamond-coated disk
High Speed Cutting
for EM2000

65E-169.03
Spare screw

65E-169.04
Special grease

65E-038.08
Screwdriver

65E-038.10
Spare 6 m cable


65E-038.11
Plastic carrying case for tooth grinding unit

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erbrich-instrumente GmbH